

Micro Fracture Chondral Pick Set

GatewayMed offers a complete set of high quality, German manufactured Micro Fracture Chondral Picks for treatment of cartilage defects on large and medium joints.

All instruments feature an ergonomic carbon fiber handle with integrated impact protection for use with the specially designed mallet.

5 different tip angles give the surgeon full flexibility in Micro fracture approaches.

Special hardened stainless steel ensures a long life-time of each Pick.



225-600-200
Chondral Pick Straight



225-600-225
Chondral Pick 25° angled tip



225-600-240
Chondral Pick 40° angled tip




225-600-260
Chondral Pick 60° angled tip



225-600-290
Chondral Pick
90° angled tip



Special hardened steel
for reliable working



Integrated Impact protection



225-600-300
Mallet for Chondral Picks

