



NEXT GENERATION OF TOURNIQUETS AND CUFFS




THE NEW 800 TOURNIQUET SYSTEM

- ENSURES MAXIMUM SAFETY
- INTUITIVE & EASY TO USE
- COMPACT & LIGHTWEIGHT



COMPREHENSIVE RANGE OF DISPOSABLE AND REUSABLE CUFFS



gatewaymed

# TOUR

Introducing the next generation in Tourniquet safety and user friendliness -

## THE 800 SERIES

Gateway Medical offers a complete range of single and double tourniquets and cuffs. With over 25 years of experience, our tourniquets are proven to be effective.



- Maximum Safety: Conforms to 3rd Edition of IEC 60601-1
- Touch Screen - No Buttons
- Easy to Use
- Simple Maintenance
- Computerized Memory
- Modern Design
- Intuitive User Interface
- Small and Lightweight
- Easy to Clean
- Possible to Customize Settings
- Printer Option
- Self Test During Set Up
- Four Models to Suit All Needs

# RNIQUETS

## TOURNIQUET MODELS

	800-20 SINGLE	800-21 SINGLE WITHOUT BATTERY	800-40 DOUBLE	800-41 DOUBLE WITHOUT BATTERY
	H x W x D	H x W x D	H x W x D	H x W x D
Dimensions (mm)	195 x 200 x 145	195 x 200 x 145	195 x 200 x 145	195 x 200 x 145
Class & Type	Class 1 Type B	Class 1 Type B	Class 1 Type B	Class 1 Type B
Line Power	100-240 VAC	100-240 VAC	100-240 VAC	100-240 VAC
Port/Cuff Connector	CPC	CPC	2 x CPC	2 x CPC
Cuff Pressure	50-550 mmHg	50-550 mmHg	50-550 mmHg	50-550 mmHg
Elapsed Time	1-16 hours	1-16 hours	1-16 hours	1-16 hours
Regulation	5 mmHg	5 mmHg	5 mmHg	5 mmHg
Weight	5.5 lbs	5.5 lbs	5.5 lbs	5.5 lbs
Battery Capacity Safe Mode	9 hours	90 minutes	9 hours	90 minutes
Battery Capacity Regulating	4 hours		4 hours	
Inflate Time	5 seconds	5 seconds	5 seconds	5 seconds
Printer Interface	USB	USB	USB	USB

## TOURNIQUET STAND



- Mobile Stand with Basket
- 5 Leg, 22" Diameter Base, for Maximum Stability
- 3" Scuff Resistant Castors for Ease of Mobility
- Safe and Convenient Storage Basket
- Adjustable Pole for Ease of Use

# CUFFS

## DISPOSABLE CUFFS



- Supplied Sterile
- Single and Double Options Available
- Latex Free
- Wide Range of Sizes
- Color Coded
- Protection Sleeve Included
- Will Fit Any Manufacturer's Tourniquet Machine
- Easy & Quick to Apply
- Highly Comfortable

### SINGLE DISPOSABLE CUFFS

Model No.	Length (cm)	Width (cm)	Length (in)	Width (in)
15-100	20	5	8	2
15-110	30	7.5	12	3
15-120	38	10	15	4
15-130	45	10	18	4
15-140	61	10	24	4
15-150	76	10	30	4
15-160	86	10	34	4
15-170	96	10	38	4
15-180	110	10	43	4

### SINGLE DISPOSABLE CUFFS CONICAL

75-140	63	9	25	3.5
75-150	75	9	30	3.5
75-160	88	11	35	4.5
75-170	110	13	43	5

### DOUBLE DISPOSABLE CUFFS

25-100	30	5	10	2
25-110	35	10	14	4
25-120	51	15	20	6
25-140	66	15	26	6
25-160	81	15	32	6

## REUSABLE CUFFS



- Machine Washable or Steam Autoclaved
- Comprehensive Range of Single and Double Cuffs
- Extra Long Silicon Bladder for Correct Application
- Velcro Tape & Strap with Plastic Clamp for Increased Safety
- Color Coded Cuffs and Hoses
- Latex Free

### SINGLE MACHINE WASHABLE & AUTOCLAVABLE CUFFS

Model No.	Item	Length (cm)	Width (cm)	Length (in)	Width (in)
35-100	Baby	39	4	15.4	1.6
35-110	Child	47	6	18.5	2.4
35-115	Arm Special	51	10	20.1	3.9
35-120	Arm	53	8	20.9	3.1
35-130	Arm Conical	65	12	25.6	4.7
35-140	Thigh Conical	78	12	30.7	4.7
35-145	Thigh Conical Special	86	14	33.9	5.5
35-150	Thigh Long Conical	93	13	36.6	5.1
35-160	Thigh Extra Long Conical	102	13	40.2	5.1
35-170	Thigh Jumbo Conical	126	13	49.6	5.1

### DOUBLE MACHINE WASHABLE & AUTOCLAVABLE CUFFS

45-90	Baby Small	31	7	12.2	2.8
45-100	Baby	41	7	16.1	2.8
45-105	Baby Wide	46	10	18.1	3.9
45-110	Child	50	10	19.7	3.9
45-120	Arm Conical	67	17	26.4	6.7
45-130	Arm Conical	82	17	32.3	6.7
45-140	Thigh Conical	100	17	39.4	6.7
45-150	Thigh Conical	125	17	49.2	6.7